

	Пробирки PCR-02D-C	Пробирки PCR-02-C	Пробирки SSI-3225	Стрипы PCR-0208-C + крышки PCR-02CP-C	Стрипы PCR-0208-C + крышки PCR-2CP-RT-C	Стрипы SSI-3240	Стрипы SSI-3247	Стрипы SSI-3248
АНК-16/32	+					+		
RotorGene 2000/3000		+	+					
RotorGene 6000	+	+	+					
iQ/iQ5	+		+	+	+	+	+	
CFX96	±		±	±	±	±	±	+
ABI Prism 7000/7300	+		+	+	+	+	+	
ABI Prism 7500	+		+	+	+	+	+	